

Call and reset button with nurse call socket

Art.: ERTH0013-1



Call and reset button with nurse call socket

Description

What is a Call and reset button with nurse call socket?

Call and reset button with nurse call socket is a wired bedside nurse call interface that connects a patient's call cord to the room terminal. This version includes a red call button for initiating a nurse call and a green reset/presence button used by staff to register their presence when attending the patient and to reset active calls. The unit ensures reliable, low-latency communication, enabling quick and confident patient-caregiver interaction.

Features and Benefits

- **Wired Connection Through the Whole System:** Ensures fast, stable, and interference-free communication between bedside units and the room terminal.
 - **Red Call Button with LED:** Allows patients to quickly and easily trigger a nurse call.
 - **Green Reset/Presence Button with LED:** Enables staff to confirm their presence at the bedside and reset active calls directly at the point of care.
- **Call Cord Disconnection Detection (Pull Alarm):** Automatically triggers an alert if the nurse call cord is unplugged, ensuring safety and immediate attention.
- **Compact, Hygienic Design:** Plastic housing with a waterproof, antibacterial membrane keypad supports easy cleaning and infection control compliance.

- **Compliance & Quality Assurance:** Fully compliant with DIN VDE 0834 standards and CE certified, ensuring top safety and performance standards.

Best Settings

- ICU
- General Wards
- Emergency Departments
- Post-Operative Care
- Assisted Living Environments
- Pediatric & Neonatal Units

Technical specification

Power supply

- Power-In: 24V DC
- Power Consumption: 5mA

Dimensions

- W:85mm / H:85mm / D:26mm

Weight (Net)

- 25g

Material

- Plastic housing
- Waterproof, antibacterial membrane keypad

Operating Conditions

- Operating Temperature: +5°C to +55°C
- Operating Humidity: Up to 85%, non-condensing

Other

- Wired connection to the call cord via RJ45 port
- Indicator on Membrane: Positioning indicator (LED) for night visibility of the device / optical call display

- Mounting: Wall-mounted (flush or surface installation)

Caution

- Without holder and frame - to be bought separately!

Compliant and tested with CE certification requirements.

Compliant and tested with DIN VDE 0834 requirements.

*Disclaimer: All product specifications, including dimensions, weight, features, and performance, are subject to change. To ensure accuracy, please verify all details before purchase.