

Wired call button for gentle call activation

Art.: ERTH0017



Wired call button for gentle call activation

Description

What is a Wired call button for gentle call activation?

The wired call button for gentle call activation is a call button for gentle call activation for individuals with visual impairments and/or physical disabilities that require a larger target area.

Features

- **Large Activation Surface:** A 115 mm activation area.
- **Reliable Sensory Feedback:** Provides both tactile (a click) and auditory (sound) feedback.

Benefits

- **Enhanced Accessibility:** easy activation of a call with a large activation area.
- **Clear Feedback:** provides clarity of action.

Best Settings

- General Wards.
- ICUs.
- Post-Operative Care.

- Pediatric & Neonatal Units.

Technical specification

Dimensions

- L:16.2 cm / W:16.2 cm / H: 3.5 cm*

Weight (Net)

- 275 g*

Material

- Plastic housing

Operating Conditions

- Operating Temperature: +5°C to +40°C
- Operating Humidity: Up to 85%, non-condensing

Other

- Activation force: 150 g
- Activation area: 115 mm
- Auditory feedback: Click on activation
- Tactile feedback: Yes
- Visual feedback: No
- Mounting: Quick Ready Mounting Plate for Assistive Technology
- Connector type: RJ45
- Cable length: 2.5 m

Compliant and tested with CE certification requirements.

Compliant and tested with DIN VDE 0834 requirements.

*Disclaimer: Product dimensions and weight are subject to change. Please verify details before purchase.