

Wireless motion detector with call-reset button

Art.: ERTW0067-1



Description:

Wireless motion detector with red call button and green presence button. Usually the device is placed under the bed and when motion is detected, the device sends SOS to the NurseCare interface / Room terminal.

Technical specification:

- Wireless call/reset button (Red and Green)
- LED indicator: red for call button, green for reset button
- Waterproof membrane keypad
- Frequency band: 2.4GHz
- Voltage: battery 3V

- Battery type: CR2450
- Dimensions (W:85mm / H: 85mm / D:10mm)
- Weight: 25g
- Compliant with CE certification requirements
- Compliant with DIN VDE 0834 requirements

PIR Sensor Characteristic:

- Detection distance: 5 meters
- Detection range: 94 degrees x 82 degrees
- 64 detection zones